



(01)00850043803003(10)20613

<b>Catalog Number</b>	PSPD-1000
<b>Description</b>	Planet-Safe <sup>®</sup> Petri Dishes, 100mm x 15mm, Sterile, 375/case
<b>Manufacture Date</b>	June 10, 2022
<b>Lot Number</b>	20613
<b>Sterilization Lot Number/Date</b>	39230/December 3, 2022
<b>Country of Origin</b>	United States of America

**Quality System Compliance**

Manufactured under FDA cGMP guidelines as presented on the date of manufacture listed above. This product is manufactured in a controlled environment, fully automated and without human intervention of parts contained within the main sterility barrier. All products are inspected during production according to applicable specifications and have passed quality inspection procedures prior to final approval for release.

**Product Specifications**

This product is a Class 1 Device as defined by the FDA.

**Material:** Resin used to manufacture this product does not contain latex, BPA, DEHP or any Phthalates.

**Sterility:** Product has been sterilized and dosimetrically released per the requirements of ISO 11137 and meets a minimum Sterility Assurance Level (SAL) of  $10^{-3}$ .

**BSE/TSE:** Product does not contain, make use of, or involve any raw materials of animal origin, including animal proteins. No products of animal origin are stored at the manufacturing facility.

**Non-Pyrogenic:** This product has been manufactured in a contained and controlled manufacturing cell designed to maintain very low bioburden conditions. Therefore, all surfaces within the primary sterility barrier should be free of pyrogenic contaminant.

**DNase/RNase Free:** Resin in virgin form is FDA compliant and is molded at high temperatures without the use of additional release agents or additives allowing products to remain DNase/RNase free.

*Products manufactured/irradiated under the above lot numbers have met all quality standards at the time of release.*

**Gregory D. Kwak**  
Manager - Quality/Compliance  
Diversified Biotech, Inc.

Creation Date: 2022-01-15